

**Department of Physics and Astronomy
University of Iowa
29:137 Astronomical Laboratory
Fall 2009
Lab 4: Stellar Spectroscopy**

1 Introduction

Throughout your astronomy education, you have read about stellar spectra. Now is your chance to actually take some and compare them with the spectra presented in textbooks. We will mount the DSS-7 spectrograph on a Celestron 8 and take spectra of a few bright stars.

1.1 Reading Material

Once again, the manual for the DSS-7 spectrograph (on the course web page under “Textbooks”) will be indispensable for this project. Among other things, it describes how to get a star on the slit so we can do spectroscopy. Once the data are taken, you will need to refer to your introductory astronomy book for information on stellar spectral classification.

2 Equipment to be used

The following pieces of equipment will be used in this lab. We have equipment for two complete setups as described below. The intention is that the spectrograph has been left intact since its last calibration. You should consult with the Working Team that carried out this calibration for the data that will permit you to measure wavelengths.

- A “nosepiece” that connects the DSS-7 spectrograph to an adapter ring which fits a Celestron 8 telescope. The nosepiece screws into the slit aperture of the DSS-7.
- The ST-402XME CCD camera
- Power unit and USB cable for the ST-402XME

- The DSS-7 Deep Space Spectrograph
- A 7 inch phoneline connecting the spectrograph and the camera
- Adapter ring to connect the DSS-7 and ST-402XME
- A set of Allen wrenches, including little ones
- A laptop computer which has the program CCDOPS installed, and had drivers activated. This computer will control the spectrograph and take data.

3 System Setup

The spectrograph was assembled in the previous class period. Assume that this was done right and that the wavelength calibration is still valid. Carry all the equipment to the roof, and set up the Celestron 8 telescopes. This should be done while there is still daylight. The whole system should be assembled and turned on while there is still light to work by. We can then wait until it gets dark to take spectra.

Follow these steps.

1. Screw the nosepiece into the DSS-7, then connect the nosepiece to the adapter ring for the Celestron 8. Attach the adapter ring to the C-8. There may be real problems with the weight distribution on the Celestron 8. The combined camera and spectrograph are sufficiently heavy to possibly cause problems in pointing and tracking. We'll see.
2. When installing the spectrograph, try and place it so you have free access to the finder telescope and the focus knob. Both will be necessary.
3. Connect the phone line between the DSS-7 and ST-402XME, and turn on the DSS-7 with its own switch.
4. Connect the USB cable from the ST-402XME to the control computer.
5. Connect the power unit for the ST-402XME.
6. Bring up CCDOPS and establish a link between the computer and the spectrograph. Listen for the contented clicking noises of the ST-402XME.

4 Putting the Star on the Slit and Taking Spectra

We will start by practicing with Vega (Alpha Lyrae). Taking stellar spectra consists of closely spaced operations of manoeuvring the star onto the slit, then quickly taking a spectra (“Grab Spectrum”) before the star drifts out of the slit. This will probably be a difficult and exasperating operation in view of the coarse position controls on the Celestron 8 telescopes. Be prepared to take a lot of time doing this. Be patient. Curses, blood-curdling oaths, and weeping outbursts of self-pity will not be tolerated in the laboratory. Follow this procedure.

1. Go to the section of the DSS-7 manual entitled “Operation at the Telescope”. Follow the instructions, as summarized below.
2. A crucial task is to place a marker in the field of view that indicates the slit location. This will permit you to guide the star onto the slit. First, be sure the box “DSS-7 Enabled” is checked in the DSS-7 menu box. Then choose VIEW SLIT, then click POSITION MARKER. This will place a white rectangle on the screen. Change the size and reposition it so it fits right over the narrowest slit.
3. You are now ready to position your star on the slit, in preparation for measuring its spectrum. Click on the DSS-7 dialog box, and select “POSITION MODE”. Choose a short integration time and take an image. As described in the manual, this is like an image frame, except for the little white box which indicates the position of the slit. While observing the image of the star, move the RA and DEC controls to get the star onto the white box¹.
4. When you have the star on the slit, bring up the dialog box, choose “GRAB SPECTRA”, and take a spectrum. Start with an integration time of 10 seconds. Check the resultant spectrum to see if it is usable (i.e. not saturated, but with an adequate signal-to-noise ratio). Take a number of spectra with different exposure times. If the ST-402XME

¹Curses, blood-curdling oaths, and weeping outbursts of self-pity will not be tolerated in the laboratory.

is satisfactorily aligned with the DSS-7 so that the spectra are parallel to the horizontal axis, you could turn on the “vertical binning” at 4 pixels to increase the signal-to-noise ratio.

5. Save your spectra to a folder for later processing. Choose file names that tell you what star and exposure time was used. If further notes are necessary, one member of the Working Team can act as scribe, and share notes with other Team members later.
6. In order to permit analysis with the SBIG spectroscopy analysis program, you should also save a copy of the spectrum that is “cropped”. There is a utility in CCDOPS that cuts out a 20 pixel-high ribbon that can be read by the analysis program DSS7 (referred to as SCP in the manual). In CCDOPS, go to the UTILITY menu, then select CROP, then select SPECTROSCOPY. The program will pop a wide, white band on your set of spectra. Position it so it fits right on the spectrum you want to save (usually the one with the narrowest slit, and thus the highest spectral resolution). When you have positioned the CROP mask, go back to the UTILITY menu, select CROP, and under the choices, pick CROP IMAGE. Pick a name for the file that identifies it as a cropped spectrum, then save to your folder of spectra for the night.

5 Observations

If the observations work properly, take spectra of the following stars. The instructor may alter the selection of stars, depending on sky conditions, time of night, and degree to which the equipment is functioning.

1. Vega (Alpha Lyrae)
2. Arcturus (Alpha Bootes). This is the bright star in the western sky.
3. Altair (Alpha Aquilae)
4. η Cassiopeiae (Eta Cassiopeiae). This is a G3V star at a distance of 5.94 parsecs. As you will recognize from its spectral type, it is nearly the same type of star as the Sun. It is easy to find in the sky, but it is

several magnitudes fainter than the other stars we have observed, so a longer integration time will be necessary.

5. The final measurement is taking of a calibration spectrum, so we can determine the wavelength scale for these observations. Pick a street light for observations. Be sure and choose a mercury or sodium vapor light rather than an incandescent one. The lights that appear bluish-green are mercury, and those that are yellow are sodium. Pick a sodium light that has a purer yellow light, rather than the pale light. The former class are low pressure sodium lights, and they have a less pressure-broadened spectrum. Save the calibration spectra files in both fits and cropped format.

6 Analysis and Results

The preferred approach to the analysis is with the SBIG spectroscopy program DSS7 (as it appears in the SBIG directory on the computers in the lab), referred to in the manual as the Spectral Calibration Program (SCP). A complete description of the use of this program appears in the section of the manual entitled “Use of the SBIG DSS-7 Spectral Calibration Program (SCP) software”. Follow the instructions in the manual. One point to make at the outset is that the program is designed around reading in both a calibration spectrum (our street lights) as well as the spectrum to be studied (the stellar and planetary spectra). You should choose your best calibration spectrum before starting the session (or you can blunder through, learn from your mistakes, and do it right the second time).

Analysis of the spectra can also be done with Maxim DL. A slice can be made along the horizontal axis of the image corresponding to a spectrum, which gives a plot of intensity versus pixel position. This can be converted to wavelength with the calibration curve developed in Lab 3. there are five slits on the DSS-7, giving different spectral resolutions. Use the narrowest slit (best spectral resolution) if possible, but move up to coarser spectral resolution if you need more signal.

The best way to carry out the analysis tasks described below is with the SCP. **However**, I don’t want you to just accept the wavelength calibration as something divinely-revealed. You can and should analyse the calibration

spectra from scratch, using Maxim DL, and verify the wavelength calibration that SCP gives you. I will look for this in your report. The SCP has options for generating and printing out line drawings of the spectra.

1. Make plots of the spectra for each of the stars for which you have results. Include these with the worksheets that you hand in.
2. For each star, measure the wavelength of the spectral maximum, and use Wien's Law to estimate the photospheric temperature of the star. In some cases, this may be difficult, or you may only be able to provide an upper limit or lower limit. Clearly indicate this. **Be organized!** Set up your measurements and data in clearly labeled tables.
3. Identify absorption (or emission) lines in the spectrum of each star. Using the measured wavelengths, identify these spectral lines. This will probably involve consulting an astronomy book which contains information on stellar spectra.
4. On the basis of item # 2 and # 3, determine the spectral class of each star. Look up the "real" spectral class and note that down, too.
5. Estimate the wavelength resolution in your spectra. This may differ from one star to another.
6. Save your spectra files. After I have had a chance to examine them, I may propose additional analysis tasks that you could undertake in future lab sessions.